

# PQuery

- **PQuery** is a flask web application (prototype) built to visualise annotated multi-sample variant-level data stored in VCF files.
- It's really fast for querying small regions like multiple genes.
- User can select genomic regions and samples to query:
  - If  $\leq 40$  samples selected, PQuery will run in **sample mode** and the output table will include per-sample genotype information.
  - If  $> 40$  samples selected, PQuery will run in **cohort mode** and the output table will not include per sample genotype columns in the output table.
- It is better to query small regions rather than the whole exome which will be slow.
- PQuery is querying the hg19 POlab data.

# Home Page

PQquery About

PQquery is an application developed by the **Genetics Centre** of St George's University of London for fast and efficient interpretation of the Lymphoedema and Lipoedema Whole Exome Sequencing Data Inventory consisting of more than 250 samples.

**Samples Selection**

Select Sample None selected OR Select Cohort None selected OR Upload CSV file Browse... Info

**Genomic Region Selection**

Select Genomic Region None selected OR Enter Gene Name or Genomic Coordinates EZH2,ENG00000012,chr5:11 OR Upload CSV file Browse... Info

Submit

Dropdown menu with all samples (with embedded search box)

Dropdown menu cohorts (Lymphoedema or Lipoedema)

Free form text box. For info hover over the Info button

Update CSV file with sample names or genomic locations (one per line) – Hover over the info buttons on each field for more info.

- You have to use **only one** of the proposed ways to enter your input for samples and/or for genomic regions.
- You need to specify at least one sample and at least one valid genomic regions or the app will throw an error.
- Hit submit when you enter your input.
- Everytime you hit submit (or refresh) the page is being refresh and the input is being cleared.

# This is the var\_table page with the query result: An ultra-wide annotated variants table

Download the table as CSV button.

CSV

This is the only way you can filter the table:

1. Select the to-be-filtered column.
2. Specify your filters.
3. Then you hit Submit.

See next slide.

Filtering Pane

Select the Columns you would like to filter:

Select Columns...

Submit

Click the Inspect IGV button to check the IGV for the selected samples

About

Inspect IGV

Columns showing the Alt. Allele Count (AC), Total Allele Number (AN) and Allele Frequency (AF) for the **SELECTED** samples only. They are **dynamically** calculated.

Plain AC, AF, AN columns correspond to the whole cohort.

variant_id	chr	start	end	ref	alt	qual	filter	ac	af	an	dynamic_ac	dynamic_an	dynamic_af	func.refgene	gene.refgene	genedetail.refgene
187594	7	100403121	100403121	A	C	75.78	VQSRTTrancheSNP99.90to100.00	6	0.011	556	1	12	0.0833333333333333	splicing	EPHB4	NM_004444:exon15:c.2678+2T>G
187600	7	100405105	100405105	C	T	8123.09	PASS	5	0.008993	556	1	12	0.0833333333333333	exonic	EPHB4	.
187603	7	100410597	100410597	G	A	393949	PASS	174	0.313	556	4	12	0.3333333333333333	exonic	EPHB4	.
187604	7	100410657	100410657	A	G	867977	PASS	374	0.673	556	8	12	0.6666666666666666	intronic	EPHB4	.
187606	7	100411278	100411278	T	C	592542	PASS	374	0.673	556	8	12	0.6666666666666666	exonic	EPHB4	.
187608	7	100411504	100411504	A	AG	2784.84	PASS	3	0.005396	556	1	12	0.0833333333333333	intronic	EPHB4	.
187615	7	100414788	100414788	A	C	90.15	VQSRTTrancheSNP99.90to100.00	7	0.013	556	1	12	0.0833333333333333	intronic	EPHB4	.
187618	7	100416250	100416250	A	G	781325	PASS	374	0.673	556	8	12	0.6666666666666666	exonic	EPHB4	.
187621	7	100420155	100420155	A	G	67765.6	PASS	47	0.085	556	1	12	0.0833333333333333	exonic	EPHB4	.
187626	7	100421781	100421781	C	T	55222.4	PASS	48	0.086	556	1	12	0.0833333333333333	intronic	EPHB4	.
187628	7	100424577	100424577	G	GC	23806.8	PASS	78	0.148	526	2	12	0.1666666666666666	intronic	EPHB4	.
187630	7	100424586	100424586	CC	C	113500	PASS	160	0.295	540	4	12	0.3333333333333333	intronic	EPHB4	.

Showing 1 to 26 of 26 entries

Previous 1 Next

Use your browser's find functionality (Ctrl+F) to search for column names

## Entries counter: BUG

This is **not** showing the right number of variants returned.  
It shows the number of variants BEFORE they are filtered for AC>0. Ignore for now.

Scroll to the very left end of the table to see the **genotypes per variant of the selected samples IF RUNNING ON THE SAMPLES mode**.

Use the normal Scrollbars to navigate through the table in both dimensions.

Paging options

## Filtering Pane

Inspect IGV

Select the Columns you would like to filter:

D102715

ref

func.refgene

ncRNA\_exonic  
ncRNA\_intronic  
exonic  
intergenic

gene.refgene

cadd\_phred From:

To:

D102715

Select...

Submit

There are 4 types of filter predefined for each column:

1. Free form text (For example for the REF/ALT columns).
2. Multi-select dropdown menu.
3. Numeric Ranges (From, to) If you want you can only specify 'from' that equals to '>=' or 'to' that equals to '<=' leaving the other one empty.
4. Sample Genotypes filter.

CSV

variant_id	chr	start	end	ref	alt	qual	filter	ac	af	an	dynamic_ac	dynamic_an	dynamic_af	func.refgene	sample	phred
187603	7	100410597	100410597	G	A	393949	PASS	174	0.313	556	3	10	0.3	exonic	EPHB4	.
187604	7	100410657	100410657	A	G	867977	PASS	374	0.673	556	7	10	0.7	intronic	EPHB4	.
187606	7	100411278	100411278	T	C	592542	PASS	374	0.673	556	7	10	0.7	exonic	EPHB4	.
187615	7	100414788	100414788	A	C	90.15	VQSRTTrancheSNP99.90to100.00	7	0.013	556	1	10	0.1	exonic	EPHB4	.
187618	7	100416250	100416250	A	G	781325	PASS	374	0.673	556	7	10	0.7	intronic	EPHB4	.
187621	7	100420155	100420155	A	G	67765.6	PASS	47	0.085	556	1	10	0.1	exonic	EPHB4	.
187626	7	100421781	100421781	C	T	55222.4	PASS	48	0.086	556	1	10	0.1	intronic	EPHB4	.

Click **Submit** to apply the filters. Whenever you click submit all the specified filtering is being cleared and YOU NEED TO RESPECIFY all of them if you want to add one more filter for example.

Use these small arrows to **sort** the columns. You can do that on **every column** on the table

Select the Samples you would like to check IGV for:

None selected

And then click Submit:

Submit

Click here to show the  
list of samples which  
can be selected for IGV.

Select the samples.

Then click Submit.

After clicking submit  
you can set-up IGV. For  
more information on  
how to use IGV  
WebApp check here:  
<https://igvteam.github.io/igv-webapp/>

